



**Surface Measurement Systems**  
World Leader in Sorption Science

# Sorption technology for pharmaceutical applications

Nachal Subramanian

Surface Measurement Systems



# Overview

- Surface Measurement Systems - Introduction
- Dynamic Vapour Sorption (DVS)
- Inverse gas chromatography (iGC)

# Surface Measurement Systems

**Surface Measurement Systems** develops and engineers techniques and instrumentation for physico-chemical characterisation of complex solids. We are the world leaders in Dynamic Vapor Sorption technology and Inverse Gas chromatography instrumentation, providing professional world-class scientific and technical support for our international customers.

Our range of characterization instruments continues to help solve difficult problems in the pharmaceuticals, biomaterials, polymers, catalysts, chemical, cosmetics, building materials and food industries, and are used by hundreds of leading laboratories and universities throughout the world.

# Overview


- Surface Measurement Systems - Introduction
- **Dynamic Vapour Sorption (DVS)**
- Inverse gas chromatography (iGC)

# Dynamic Vapour Sorption

## What is DVS?


DVS is a gravimetric technique that measures how quickly and how much of a solvent is absorbed by a sample.

- Sample is exposed to a series of step changes in relative humidity and the mass change is measured as a function of time.



# What can the DVS do for me?

1. How does my material interact with moisture or solvents and temperature in the vapour phase?
2. **Stability, Performance and Processing issues: Reversible and Irreversible effects of moisture**
3. **Create Moisture Isotherms** – i.e. Equilibrium moisture content as a function of %RH
4. **Heterogeneity?** – Identify the Heterogeneity of a sample batch
5. Homogeneity? – Identify variance within one sample
6. **Kinetics** – Moisture transport properties, how fast or slow?
7. Energy – How strongly is the moisture bound to the material, surface or bulk?
8. **Identify & Characterise Phase Transition/Changes, e.g. polymorphs, amorphous stoichiometry**
9. Hydration and Solvate Formation
10. **Drying Analysis**
11. Diffusion and Activation Energy
12. **Heat of Sorption**
13. Moisture Uptake/Content? i.e. how much moisture/vapour is taken up or release



# Comparison to traditional jar method

## Jar Method (static sorption)




- Static system - Very slow achievement of equilibrium
- Longer experiment time
- Large amounts of sample required (1 to 10 g)
- Risk of contamination or sample loss due to manual weighing

## DVS (dynamic sorption)



- Dynamic Flowing Gas System - Faster Equilibrium
- High sensitivity Ultra-Microbalance (0.1  $\mu\text{g}$  resolution)
- Allows small sample to be used (1 to 10 mg)
- No risk of contamination or sample loss
- Both sorption and desorption measured

# Typical DVS data










# Applications of DVS - Examples

## Phase transitions - Amorphous Lactose

DVS Data




Raman Data



# Applications of DVS - Examples


## Glass Transition and Crystallisation




# Pure Sucrose Images



□ 0% RH



□ 32.7% RH Glass Transition



□ 80.9% RH Crystallized



□ 92.6% RH Deliquesced

# Properties measured by DVS

**DVS is a tool for thermodynamic and kinetic studies of surface and bulk properties**

- Sorption Isotherms
- BET Specific Surface Area
- Phase Transitions
- Permeability and Diffusion
- Competitive (Multicomponent) Adsorption
- Kinetics information
- Heat of Sorption
- Tg RH determination
- Amorphous content determination
- Camera and Raman capability

# DV<sup>c</sup> Family of Products



**DVS**  
**INTRINSIC PLUS**

Water vapor only  
20-40 °C  
Small footprint



**DVS**  
**VACUUM**

10<sup>-6</sup> Torr to Atmospheric  
20-70 °C  
**Water + Organic Vapors & Gases**  
Competitive sorption  
Sample preheat to 400 °C



**DVS**  
**ADVENTURE**

Water only  
**5-85 °C**  
Temp stability +/-  
0.1°C  
Camera and Raman  
Sample preheat up to  
300 °C



**DVS**  
**RESOLUTION**

**Water + Organic Vapors & Gases**  
**5-85 °C**  
**Speed of Sound Sensor**  
Temp stability +/- 0.1°C  
Camera and Raman  
Sample preheat up to  
300°C  
Has 2 configurations:  
Standard and Advanced

**DVS**  
**ENDEAVOUR**



**5 experiments in parallel**

# Overview

- Surface Measurement Systems - Introduction
- Dynamic Vapour Sorption (DVS)
- **Inverse gas chromatography (iGC)**



# Inverse Gas Chromatography

- Gas phase sorption technique.
- Focus on physicochemical studies - kinetic information and thermodynamic quantities from sorption equilibria.
- Earlier work on catalytic materials, e.g. activated carbon, alumina and silica.
- Powerful physico-chemical characterization tool for **powders, fibres, films, particulates, semi-solids**.
- Surface Energy is the most common measured property by iGC.


\*\* Schultz, J., Lavielle, L., and Martin, C., The role of the interface in carbon fibre epoxy composites. J. of Adhesion, 1987. 23(1): p45-60.

# IGC Principles

## Analytical Gas Chromatography



## INVERSE Gas Chromatography





# iGC-SEA Introduction

- **Gas phase injection** (like Headspace) - 12 vapor reservoirs (50 ml)
- **Carrier gas** is helium/nitrogen
- 2 column position oven design: **20 to 150 °C**
- **Background Humidity Controller**
- Flame Ionization Detector (FID)
- **User Friendly** Control and Analysis **Software**



Carbon fiber   cotton   hair   granules   powder   medical metal implants



**Safety Features:**  
**Hydrogen Leak & Organic Vapor Leak Detectors**

# Properties measured by IGC-SEA

The IGC-SEA provides unique access to the following physico-chemical properties of a wide range of solid materials in a controlled humidity environment:

- Dispersive and Polar **Surface Energies**
- Heats and Entropies of Adsorption
- **Acid/Base** Interactions
- **BET** Specific Surface Area
- Phase Transitions
- Sorption **Isotherms**
- Permeability, Solubility and Diffusion
- Competitive (Multicomponent) Adsorption
- Thermodynamic Work of **Cohesion** and **Adhesion**
- Surface Energy **heterogeneity mapping**
- **Constantly extend the applications – future applications e.g. Chemisorption**

# Thank you!

Head Office:

Surface Measurement Systems, Ltd  
5 Wharfside, Rosemont Road  
London HA0 4PE, UK

Tel: +44 (0)20 8795 9400

Fax: +44 (0)20 8795 9401

Email: [info@surfacemeasurementsystems.com](mailto:info@surfacemeasurementsystems.com)

United States Office:

Surface Measurement Systems, Ltd, NA  
2125 28th Street SW, Suite 1  
Allentown PA, 18103

Tel: +1 610 798 8299

Fax: +1 610 798 0334